5.15 Radionuclide Identification

High Purity Germanium (HPGe) is a semiconductor with a bandgap of $E_g = 0.66 \text{ eV}.[20]$ At low temperatures, the density of electrons in the conduction band (and the density of holes in the valence band) in Ge is low enough that it is essentially an insulator. A high voltage across a crystal of germanium will generate an electric field in that crystal, but no current will flow.

If the crystal is then struck by a photon with sufficient energy to promote an electron from the valence band to the conduction band, the result will be an electron-hole pair. The hole and the electron will be swept in opposite directions by the electric field, and will register as a current pulse when they reach the electrodes creating the high voltage. If all of the photon energy is deposited in the crystal, then the number of electron-hole pairs generated by the incoming photon depends on the photon energy:

$$N = \frac{E_{\gamma}}{E_g}$$

So it is possible to use the magnitude of these pulses to measure the energy of the incident photons.

Warning!

The HPGe detector detects radiation by observing the small current spikes of electron-hole pairs created in the Ge crystal by the incident radiation. There is a high voltage across the crystal, which is normally non-conducting. If the crystal reaches a 'high' temperature, the thermal energy generates electronhole pairs which make the crystal conductive, the high voltage creates a high current (briefly) in the crystal, and we have to buy a very expensive new detector. Make *absolutely sure* that there has been liquid nitrogen in the detector dewar for at least six hours before you turn on the high voltage.

Procedure

- 1. Use the known sources and the chart of the nuclides (or other reference source) to make a plot of energy vs MCA channel number. From this plot, determine an appropriate calibration equation for converting channel number to energy.
- 2. Find the γ energy(ies) of the unknown sources, and use these energies to identify the unknowns. Bear in mind that many radioactive sources

decay via a long chain of reactions, each with a characteristic half-life. Depending on the age of the sample relative to its half-life, you may find that there are numerous daughter isotopes in the sample, which can complicate your analysis.